





HALF GEL COATED BODY SLEEVING



Excellent for
skin grafts,
burns, scars &
Lymphedema



- 10 Foot Rolls – 3”, 4”, 5” Flat Width
- Cut to length- will not fray
- Coated ½ side with Silipos® proprietary gel, breathable non site side
- New Fabric-Nylon/Spandex-will not stretch out
- Good for hard to reach sites
- Skin Grafts, Burns, Scar Management, Lymphedema, & Skin Care

HALF GEL COATED BODY SLEEVING

Circumferential Measurements

| Product No. | ~ Roll Size | ~ Flat Width | ~ Relaxed | ~ Fully Stretched |
|-------------|-------------|--------------|-----------|-------------------|
| 10925 | 3” x 10’ | 3” | 5.5” | 16” |
| 10935 | 4” x 10’ | 4” | 7.5” | 18” |
| 10945 | 5” x 10’ | 5” | 9.5” | 20” |

