

## NEW PRODUCT ANNOUNCEMENT

**DuraGel™ Suspension Sleeve offers a durable and skin-friendly solution.**

- Skin friendly gel works well with suction application
- DuraGel mineral oil gel moisturizes and conditions skin while providing increased durability
- Available in 3mm gel thickness
- Easy to take on and off



<i>Basic</i>		<i>Reinforced</i>	
<i>Part</i>	<i>Size</i>	<i>Part</i>	<i>Size</i>
18181	Medium	18191	Medium
18182	Medium +	18192	Medium +
18183	Large	18193	Large
18184	Large +	18194	Large +
18185	Medium + Distal +	18195	Medium + Distal +

### DuraGel™ Suspension Sleeve Sizing Guidelines

Size	Proximal Circumference		Distal Circumference	
	Relaxed	Fully Stretched	Relaxed	Fully Stretched
Medium	9" (23cm)	15" (38cm)	8" (20cm)	13" (33cm)
Medium+	14.25" (36cm)	26.75" (68cm)	9.25" (23cm)	16" (40cm)
Medium+ Distal+	14.25" (36cm)	26.75" (68cm)	9.75" (24.5cm)	17" (43cm)
Large	13.5" (34cm)	24" (61cm)	12" (31cm)	18" (46cm)
Large+	16.25" (41cm)	28" (71cm)	11.25" (29cm)	20" (51cm)

*All sizes are approximate.*