

SILIPOS®

SiloLiner™

Sizing Guidelines: Pin Liners


SiloLiner™ Product No.	Flat Width	Distal Circ. at urethane attachment	Proximal Circ. MPT
1460	3" / 7.6cm	8-12" / 20-30.5cm	8-14" / 20-35.5cm
1461	4" / 10.2cm	9.5-13.5" / 24-34cm	9.5-16" / 24-40.5cm
1462	5" / 12.7cm	12-15" / 30.5-38cm	12-20" / 30.5-51cm
1463	6" / 15.2cm	14-16" / 35.5-40.5cm	14-26" / 35.5-66cm

Sizing Guidelines: Cushion Liners

Cushion SiloLiner™ Product No.	Flat Width	Distal Circ. at 2" from Distal End	Proximal Circ. MPT
1612	3" / 7.6cm	8-14" / 20-35.5cm	8-14" / 20-35.5cm
1613	4" / 10.2cm	9.5-16" / 24-40.5cm	9.5-16" / 24-40.5cm
1614	5" / 12.7cm	12-17" / 30.5-43cm	12-20" / 30.5-51cm
1615	6" / 15.2cm	14-19" / 35.5-48cm	14-26" / 35.5-66cm

Low Activity Liner

Sizing Guidelines




Liner Size	Proximal Circumference (Measure 15" from Distal End)			
	Proximal		Distal	
	Relaxed	Fully Stretched	Relaxed	Fully Stretched
Medium	9" / 23cm	16" / 41cm	8" / 20cm	12" / 30cm
Medium +	13" / 33cm	21" / 53cm	9" / 23cm	12" / 30cm
Large	15" / 38cm	22" / 56cm	9" / 23cm	15" / 38cm
Large +	16" / 41cm	27" / 68cm	12" / 30cm	16" / 41cm
X Large	16" / 41cm	30" / 76cm	12" / 30cm	20" / 51cm

Distal Circumference (Measure 1.5" from Distal End)

DuraGel™

Sizing Guidelines



Liner Size	Proximal Circumference (Measure 15" from Distal End)			
	Proximal		Distal	
	Relaxed	Fully Stretched	Relaxed	Fully Stretched
Medium	9" / 23cm	15" / 38cm	8" / 20cm	13" / 33cm
Medium +	13" / 33cm	20" / 51cm	9" / 23cm	13" / 33cm
Large	15" / 38cm	21" / 53cm	9" / 23cm	16" / 41cm
Large +	16" / 41cm	26" / 66cm	12" / 30cm	17" / 43cm
X Large	16" / 41cm	29" / 74cm	12" / 30cm	21" / 53cm

Distal Circumference (Measure 1.5" from Distal End)

DuraGel™ by SILIPOS

Suspension Sleeve

DuraGel™ Suspension Sleeve Sizing Guidelines

Size	Proximal Circumference		Distal Circumference	
	Relaxed	Fully Stretched	Relaxed	Fully Stretched
Medium	9" (23cm)	15" (38cm)	8" (20cm)	13" (33cm)
Medium+	14.25" (36cm)	26.75" (68cm)	9.25" (23cm)	16" (40cm)
Medium+ Distal+	14.25" (36cm)	26.75" (68cm)	9.75" (24.5cm)	17" (43cm)
Large	13.5" (34cm)	24" (61cm)	12" (31cm)	18" (46cm)
Large+	16.25" (41cm)	28" (71cm)	11.25" (29cm)	20" (51cm)

All sizes are approximate.